



Features

Specifications

◎ Main Character

- Memory customized expression sequence
- Vacuum control in 10 levels
- Auto switch off after 30 minutes
- Sufficient vacuum for dual expression
- Backflow prevention
- Rechargeable battery
- Convenient and graceful design
- 3 operating mode (express , massage , memory)
- Operating time display
- Silent

Features

Specifications

Functional Characteristics

| | | | | |
|----------------|---------------------------|------------------|----------------------|---|
| Type | Dual Electric Breast Pump | Pumping Pressure | Vacuum | -20 ~ -270mmHg |
| Display | Type | Function | Mode | Expression 10levels Massage 10levels |
| | 3" Mono LCD | | Auto Shut Off | 30min |
| Pressure Level | Numeric & 10levels | Memory | Back Flow Prevention | 1 ~ 60levels |
| Operating Time | Numeric Display(min:sec) | | | |
| Battery Level | 3levels | | | |

Power

| | | | |
|-------------|---------------------------------|----------------------|--|
| Adaptor | Input : AC 100 ~ 240V (50/60Hz) | Rechargeable Battery | 3.7V Li-ion Polymer |
| Consumption | Output : DC 12V/2A | | Operating Time : 150min Charging Time : 90min |

Standard Configuration

| | | | |
|------------------|-----|---------------------|---------------------|
| Operation manual | 1ea | Funnel & Bottle Set | 1ea (Funnel : 24mm) |
| Power Adaptor | 1ea | | |

Options

| | | | |
|------------------------|-------------|--------------------------------|--|
| Carrying Case | | Funnel & Bottle Set (for Dual) | |
| Funnel | 21mm / 27mm | Cigar Cable | |
| Mobile Charging Gender | | | |

Warranty

| | |
|-----------|-------|
| Main Unit | 1year |
|-----------|-------|

Physical Characteristics

| | | | |
|-------------------|----------------------------|-------------------|-------|
| Dimension | | Weight | |
| Main Unit | 70(W) x 45(D) x 110(H)mm | Main Unit | 370g |
| Packing | 240(W) x 140(D) x 150(H)mm | Packing | 0.9Kg |
| Carton Box (10ea) | 720(W) x 260(D) x 340(H)mm | Carton Box (10ea) | 10Kg |

Environmental Conditions

| | | | |
|-----------------------|-------------------------|---------------------|--------------------------|
| Operating Temperature | 10 to 40°C (50 ~ 104°F) | Storage Temperature | -20 to 60°C (-4 ~ 140°F) |
| Operating Humidity | 5 ~ 85% non-condensing | Storage Humidity | 0 ~ 95% non-condensing |

Certificates

| | | | |
|--------------|------------------------|-----------------|-------------------|
| Certificates | KFDA | Classifications | Class II, Type BF |
| Standard | EN60601-1, EN60601-1-2 | Waterproof | IP21 |

Head Office

Bistos Co.,Ltd

7th Fl., A Bldg., Woolim Lions Valley 5 -cha, 144-3,

Sangdaewon-dong, Jungwon-gu,

Seongnam-si, Gyeonggi-do, Korea

Tel : 82-31-750-0340 Fax : 82-31-750-0344

E-mail : bistos@bistos.co.kr <http://www.bistos.co.kr>