



1. External breast compression and artificial respiratory: indicator light display, digital display device and sound prompting device;
2. Operation mode: exercise mode; examine mode
3. Operation time: can count and set the operation time
4. Sound prompting device: volume control; turning on or turning off sound prompting device
5. Print: print operation result
6. Examination of pupil response: pupil will change into normal size after accomplishing the rescue process;
7. Examination of carotid response: touch carotid pulse by hand; simulating spontaneous carotid pulse in compression process; simulating spontaneous carotid pulse after examine operation;
8. Working conditions: output power is 24V, the manostat apply to various external power;

Alternative Accessories:

1. Trauma modules of limbs