

> Digital Ultrasound Combine Therapy (EL-SONIC, ELS-100)

Electronic + Ultrasound Therapy

EL-SONIC, ELS-100 is a combine device with ultrasound therapy and low frequency stimulator that can be easily applied to patient's application.

Features

- Various Auto Program of 4 types,
- Safe start from "0" intensity,
- Maximize of stimulation by Continuous / Discontinuous, interception/ current outputs,
- Constant Current method,
- Can offer intensity and amount to stimulate of sympathetic nervous system,
(composition and output of digital system)
- Control the low frequency signal in detail by Pulse width in order to minimize of feeling skin rejection,
- Form of output so that can use Combine and low frequency, ultrasonic,
(independently)(2 persons)

Accessories



Specifications

Items	Specification
Stimulation Method	Constant Current
Maximum Output Current	Max. 20mA / Rectangle / 500ohm
Mode Range	Press / Tap / Mass / Roll
Main Power	220V / 60Hz
Time	3, 5, 7, 10, 15, 20 set
Power Consumption	50VA
Frequency	1~500HZ
Ultrasound Frequency	1MHz+/-20%
Ultrasound Intensity	3W/cm ² under
Dimension	425(W)x203(D)x94(H)(mm)
Weight	5Kg
Electronic Shock Protect Class	1 BF Type Device
Excess current protection circuit	
Zero Start	

