

GD/CPR175S ADVANCED ELECTRONIC HALF-BODY CPR TRAINING MANIKIN



GD/CPR175S is upgraded from GD/CPR170S, with electronic display of CPR operation and the main functions; simple operation and cost-effective.

Implementation standards: Complying with AHA (American Heart Association) 2005 guideline for CPR and ECC (emergency cardiovascular care)

Features:

1. Simulate standard airway open
2. Manual external breast compression with electric light display
Correct compression position, yellow light will power on;
Correct compression intensity (4-5 cm areas), green light will power on;
Wrong compression intensity (>5cm areas), right light will power on;
3. Artificial mouth to mouth respiration (inhalation)
The volume of inhalation can be judged through observing the wave of breast (volume standard $\leq 500/600-1000\text{ml} \leq$)
4. Operation frequency: the newest international standard: 100 times/minutes
5. Operation methods: training operation