

## OXYGEN GENERATOR AND OTHER SER

# KELING MEDICAL APPLIANCES



KL-ZY5L

KL-ZY5LW type has the function of atomizing

### Features

- ⊙ Full plastic shell, handsome appearance;
- ⊙ Safety valve inside, safe and reliable in application;
- ⊙ Adopting the top product SXS lithium molecular sieve of France, through pressure swing adsorption separation process, producing high purity oxygen that accord with the standards of medical oxygen.

Items	Model NO.	KL-ZY3L(W)	KL-ZY4L(W)	KL-ZY5L(W)
Major technical parameters	Mains voltage V/Hz V/H	AC220/50	AC220/50	AC220/50
	VA Power input VA	400	400	450
	Oxygen content%	≥90	≥90	≥90
	L/min Rated oxygen generating efficiency L/min	3	4	5
	ml/min Atomization efficiency ml/min	≥0.2	≥0.2	≥0.2
	Mpa Oxygen delivery pressure Mpa	0.04-0.06	0.04-0.06	0.04-0.06
	Noise level dB(A)	≤50	≤50	≤50
	Relative humidity of operating environment%	≤80	≤80	≤80
	Ambient temperature °C	5-40	5-40	5-40
	Kpa Atmospheric Alarm for Kpa	86-106	86-106	86-106
	mm Dimensions mm	300X560X575	300X560X575	300X560X575
	Kg Net weight Kg	23	25	30

Note: W is satnd for it has the function of atomization.