


YKR001



**Multiduty Lumbar Traction Bed**

Technological Parameters:

- Power supply: 220/50HZ
- Power: 360w
- Traction length: 0-250mm
- Traction force: 10-990N (adjustable)
- Total time: 99 min (adjustable)
- Continuous traction time: 9 min (adjustable)
- Intermission time: 0-9 min (adjustable)
- Swing angle: 60°
- Angular degree: +60°~-20° (adjustable)

**Multiduty Cervical And Lumbar Traction Bed**

Technological parameter:

- Power: 360W
- Electrical source: 220/50Hz
- Traction length: 0-250mm
- Traction force: 10-990N (adjustable)
- Total time: 0-99min (adjustable)
- Continuous traction time: 0-9min (adjustable)
- Intermission time: 0-9min (adjustable)
- Swing angle: 60 degree
- Angular degree: 60 degree20 degree (adjustable)



YKR002



YKR003



**Multiduty Lumbar Traction Bed**

Technological Parameters:

- Power supply: 220/50HZ
- Power: 360w
- Traction length: 0-250mm
- Traction force: 10-990N (adjustable)
- Total time: 0-99 min (adjustable)
- Continuous traction time: 0-9 min (adjustable)
- Intermission time: 0-9 min (adjustable)