

Add: 20th Building , 3399 Kangxin Road, Pudong new Dist, Shanghai ,China
 Tel: +86-21-67691499 [web: www.pinruidental.com](http://www.pinruidental.com)



High Speed air turbine hand piece tie-in	2 sets
Low speed air motor hand piece tie-in	1 set
Three way syringe (cold & hot)	1 set each
Automatic cup filler	1 set
Saliva weak suction	1 set each
Dental operation light	1 set
Rotatable tempered glass spittoon	1 set
Purified water-supplied system	1 set
Dental stool	1 set
Inductive air-locked rotary arm system	1 set
Electrical & gassous integrated foot controller	1 set
Led X-ray film viewer	1 set
Comfortable backrest & seat	1 set
Adjustable headrest	1 set

Specifications

Voltage	110~240V, 50~60Hz
Air pressure	0.55MPa
Water pressure	0.2MPa ~ 0.4Mpa