

Gynecological Examination & Delivery Table

SG-510 : Full function is controlled by elec. hydraulic system, and has memory elevation and memory back section up/down function and moving range are adjustable.

SG-520 : Full function is controlled by elec. hydraulic system, and has memory elevation function and moving range are adjustable.

SG-530 : Full function is controlled by elec. hydraulic system, without memory elevation function.

SG-540 : Elevation and back section are controlled by elec. hydraulic system, and the seat section is immovable.



FEATURES :

1. Fully automatic electric hydraulic control system.
2. The memory device is built-in the elevation function. It can be elevated while the memory up/down point automatically. (SG-520 only)
The memory device is built-in the elevation and back section function. It can be moved while the memory up/down point automatically. (SG-510 only)
3. Of course, when special position is needed, you can touch the foot switch to adjust the table back section up/down; table top up/down; and seat section up/down. (SG-510 & SG-520 & SG-530 only)
4. The toilet paper roll under the seat. When the paper is properly set, you can get it continuously without moving the back or seat mattress.
5. The surface of the table is covered with UV resistant plastic sheet, it is wear-resistant, pollution-proof and easy for maintenance.
6. Water proof and 24V control foot switch.
7. The back section can be moved to -15°, and the seat section can be moved +20° (SG-540 not available). In emergency case, it'll make the patient's blood flow backward. In order to prevent the patient from suffering shock.
8. CAN BE OPTIONED THE AUXILIARY TABLE FOR DELIVERY AND OPERATION USE.
9. All the accessories are made by Stainless Steel.

SPECIFICATION :

- * TABLE TOP : 550mm(W) x 1200mm(L)
550mm(W) x 1750mm(L) With auxiliary table
- * ELEVATION : 600mm to 850mm (Without mattress)
(maximum high/low can be followed by order requirement)
- * BACK / BREAK : +60° to -15°
- * WAIST / BREAK : +20° to 0° (SG-540 not available)
- * POWER SUPPLY : 110V or 220V, 240V 50/60Hz.
(can be followed by order requirement)

STANDARD ACCESSORIES :

- | | |
|-----------------------------|--------|
| 1. SCREEN FRAME | 1 PC. |
| 2. ARM BOARD | 1 PC. |
| 3. LEG CRUTCHES WITH CLAMPS | 2 PCS. |
| 4. FORCEPS CUPS | 1 PCS. |
| 5. FOOT SWITCH | 1 SET. |
| 6. POWER CORD | 1 PC. |
| 7. CONTAMINATION BASIN | 1 PC. |
| 8. HIP PAD (WATER PROOF) | 1 PC. |

OPTIONAL ACCESSORIES :

- | | |
|---|--------|
| A. WAIST SIDE HANDLE WITH CLAMP | 2 PCS. |
| B. DELIVERY COLLECTION PLATE | 1 PC. |
| C. CONTAMINATION BUCKET, TUBE, & FUNNEL | 1 SET. |
| D. AUXILIARY TABLE WITH MATTRESS | 1 PC. |
| E. HEAVY TYPE LEG CRUTCHES WITH CLAMPS | 2 PCS. |

STURDY INDUSTRIAL CO., LTD.
ISO-9001 & ISO-13485 CERTIFICATED FIRM



GYNECOLOGICAL Obstetric/Operation Table Examination Table



** STURDY reserves the right to improve the equipment shown, without prior notice.**

