


sonora

ultrasoundtherapy

models **sonora**
sonora combi

ultrasound
code SNR002
ultrasound+ electrotherapy
code SNR003



sonora are a programmable driver units for US transducers with working frequencies from 1 up to 3 MHz.

sonora has realized in the followings two versions: One with automatic tuning to use ultrasound transducers in the range 1 Mhz- 3 MHz. The second one **sonora combi**, with automatic tuning to use ultrasound transducers in the

range 1 MHz – 3 MHz together with electrotherapeutic current.

Through **sonora** it is possible to obtain delivery of continuous or modulated ultrasound with programmable intensity and delivery time to fit the thermo therapeutic deep treatment for different pathologic conditions.

sonora transducer handpieces can be under water used.





dimensions in mm
weight 3,5Kgs

+ what's news

- > Wide range of preset protocols
- > Customized programs
- > Multi frequency handpieces
- > Automatic recognition of the connected accessory
- > Continuous or pulsated delivery
- > Associated with electrotherapy*

sonora combi



electrotherapy

ultrasound

Main features

Output frequency (MHz)	1.0 to 3.0
Max intensity (W/cm ²)	3.0
Max output direct current (mA)*	30
Max output pulse current (mA)*	70
Max output voltage (Vdc)*	100
Power supply (Vac/Hz)	90-240/50-60
Electrical input power (VA)	70

*only with combined US and ET version (COMB)

physioled