

# ST-12 Multi-Function Deluxe Hydraulic Operation Table



1

2



3



4



5




- 1 Sitting position
- 2 Combined with ST-260 for Orthopaedic surgery positions
- 3 Kidney surgery position
- 4 Trendelenburg position
- 5 Tilted position
- 6 Head/leg plate interchangeable



6

# ST-12 Multi-Function Deluxe Hydraulic Operation Table



## Features

1. ST-12 increase flexibility to meet different surgical requirements.
2. Smooth mobile design is suitable for emergency room operations.
3. Covered moving joints offer added safety.
4. Gas spring controlled back and foot plates provide smooth manipulation.
5. Convenient and unrestricted access for C-Arm X-ray scanning.
6. ST-12 is designed to work harmoniously with ST-260 (orthopaedic extension tractor) during orthopaedic operations.
7. Height adjustable for surgeon's access and comfort (At low heights, it can accommodate seated clinicians and such as ophthalmologists and etc.).
8. Plate's durability maintains shape even after prolonged exposure to X-ray.








## Specification:

Length:	1975 mm
Width:	500 mm
Elevation:	700~1000 mm (hydraulic foot-pedal)
Trendelenburg:	$\pm 22^\circ$ (mechanical)
Tilt:	$\pm 22^\circ$ (mechanical)
Back section:	$+65^\circ \sim -30^\circ$ (gas spring)
Head section:	$+90^\circ \sim -90^\circ$ (mechanical)
Leg section:	$+15^\circ \sim -100^\circ$ (gas spring)
Brake:	Mechanical by foot-pedal

## Standard Accessories:

1. Screen frame	1 set
2. Head section	1 set
3. Arm board	1 pair
4. Leg section	1 set
5. Mattress	1 set

Sturdy reserves the right to improve the equipment shown without prior notice.



**STURDY**  
INDUSTRIAL CO., LTD.

Email: [foreign@sturdy.com.tw](mailto:foreign@sturdy.com.tw)

Website: <http://www.sturdy.com.tw>

ISO 13485  
CERTIFICATED FIRM