


# YD01-5E(LED)



## COLD LIGHT OPERATING LAMP

### Main Feature

- Excellent cold light: adopting LED as light source, without temperature rising.
- Excellent daylight quality and CRI meet the illumination needs of the doctors.
- homogenous light spot, high-definition.
- Long service life: LED life is up to 30000 hours.
- Energy saving and environment-friendly.
- Flexible and stable balance arm suspension



### Long Life Lighting

LED light source life of up to 30,000 hours or more, with an average life of an ordinary light shadowless lamp halogen 30 times, and each LED can be individually replaced, to ensure maximum use of economy.



### Technical Data

Illuminance	$\geq 100,000 \text{ Lux/}$
Colour temperature	$4800 \pm 300 \text{ K}$
Colour reduction index(Ra)	92
Total irradiance	$310 \text{ W/m}^2 / 420 \text{ W/m}^2$
Power supply voltage	A.C. $220 \text{ V} \pm 10\%$ , $50 \text{ Hz} \pm 1 \text{ Hz}$
Size of light field	160mm ~ 220mm
Illuminance depth	500mm
Lamp head diameter	500mm
LED bulb	3.2V/1W
Emergency time	3000hours

### SHANGHAI HUIFENG MEDICAL INSTRUMENT CO.,LTD.

Add No. 7753, Chuannanfeng Highway, Dattuan Town, Pudong New District, Shanghai  
Zip: 201311

Sales department

Tel: 0086-21-58082059, 58084992, 58085105, 58086595

Fax: 0086-21-58085105

E-mail: lisa@huifeng-co.com

lisa5808@163.com

Http://www.huifeng-co.com