

Die LuxaScope Sonus Familie.

- Bruststück aus Edelstahl für optimalen Klang
- Hochleistungsglasfasermembrane mit Kälteschutzring
- ergonomisch gebogener Ohrbügel aus Edelstahl mit extra verstärkter Feder für eine längere Haltbarkeit und optimalen Tragekomfort
- verschraubte weiche Ohrölvnen passen sich optimal dem Gehörgang an
- Y-Schlauch mit 3-Wege Verbindungsstück zur optimalen Schallweiterleitung in den Gehörgang
- latexfrei
- keine Verwendung von Phthalaten (Weichmachern)
- Ersatzteilset (je ein Paar harte und weiche Ohrölvnen, Membrane und Membranring) und Namensschild

The LuxaScope Sonus Family.

- stainless steel chest piece ensures excellent acoustic quality
- high performance glass fibre diaphragm with a non chill ring
- ergonomically shaped stainless steel head piece with reinforced spring for extended shelf life and optimal wearer comfort
- screwable soft ear tips adjust to the shape of the ear effectively
- Y-tube with 3-way sound transmission connector guarantees a good resonance through optimal sound conduction
- latex-free
- no use of phthalates (DEHP free)
- spare parts set (each one pair of soft and hard ear tips, diaphragm and diaphragm ring) and name plate



LuxaScope Sonus Adult

- Doppelkopf-Stethoskop für Erwachsene
- Durchmesser Membrane 44 mm | Länge 75 cm

G1.211.114 | schwarz
G1.211.214 | navy blau
G1.211.314 | grau



LuxaScope Sonus Adult

- twin-head stethoscope for adults
- diaphragm diameter 44 mm | length 75 cm

G1.211.114 | black
G1.211.214 | navy blue
G1.211.314 | grey